

Product Datasheet

Rabbit ADIPOR1 ELISA Kit (orb1292600)

Description Rabbit ADIPOR1 (Adiponectin Receptor 1)
ELISA Kit

Reactivity Rabbit

Range 0.313-20ng/ml

Target ADIPOR1


Storage 2-8 °C for 6 months

Note For research use only

Uniprot ID **B0JDQ3**

Sensitivity 0.188ng/ml

Expiration Date Please enquire.



Rabbit
ADIPOR1(Adiponectin
Receptor 1) E...