


## Product Datasheet

### Rat Spr ELISA Kit (orb1292637)

|                        |                                           |
|------------------------|-------------------------------------------|
| <b>Description</b>     | Rat Spr (Sepiapterin reductase) ELISA Kit |
| <b>Reactivity</b>      | Rat                                       |
| <b>Range</b>           | 15.625-1000pg/ml                          |
| <b>Target</b>          | Spr                                       |
| <b>Storage</b>         | 2-8 °C for 6 months                       |
| <b>Note</b>            | For research use only                     |
| <b>Uniprot ID</b>      | <b>P18297</b>                             |
| <b>Sensitivity</b>     | 9.375pg/ml                                |
| <b>Expiration Date</b> | Please enquire.                           |



Rat  
Spr(Sepiapterin  
reductase)  
ELISA Kit