


Product Datasheet

Rat Hnrnp1 ELISA Kit (orb1292638)

Description	Rat Hnrnp1 (Heterogeneous nuclear ribonucleoprotein L) ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	Hnrnp1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	F1LQ48
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat
Hnrnp1(Heterogeneous
nuclear ribonuc...