


Product Datasheet

Rat Bcl2l1 ELISA Kit (orb1292646)

Description	Rat Bcl2l1 (Bcl-2-like protein 1)ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	Bcl2l1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P53563
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat
Bcl2l1(Bcl-
2-like
protein 1)
ELISA K...