


Product Datasheet

Rat Cpn1 ELISA Kit (orb1292664)

Description	Rat Cpn1 (Carboxypeptidase N catalytic chain)ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	Cpn1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9EQV8
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat
Cpn1(Carboxypeptidase
N catalytic ch...