


Product Datasheet

Rat RTEL1 ELISA Kit (orb1292670)

Description	Rat RTEL1 (Regulator of telomere elongation helicase 1) ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/ml
Target	RTEL1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q5RJZ1
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Rat
RTEL1(Regulator
of telomere
elongati...