


Product Datasheet

Monkey CXCL10 ELISA Kit (orb1292702)

Description	Monkey CXCL10 (C-X-C motif chemokine) ELISA Kit
Reactivity	Monkey
Range	7.813-500pg/ml
Target	CXCL10
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	A0A0D9QX95
Sensitivity	4.688pg/ml
Expiration Date	Please enquire.



Monkey
CXCL10(C-
X-C motif
chemokine)
ELI...