


Product Datasheet

Mouse Wnt10b ELISA Kit (orb1292722)

Description	Mouse Wnt10b (Protein Wnt-10b) ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/ml
Target	Wnt10b
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P48614
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Mouse
Wnt10b
(Protein
Wnt-10b)
ELISA Kit