


## Product Datasheet

### Mouse Ripk3 ELISA Kit (orb1292726)

|                        |                                                                                |
|------------------------|--------------------------------------------------------------------------------|
| <b>Description</b>     | Mouse Ripk3 (Receptor-interacting serine/threonine-protein kinase 3) ELISA Kit |
| <b>Reactivity</b>      | Mouse                                                                          |
| <b>Range</b>           | 0.156-10ng/ml                                                                  |
| <b>Target</b>          | Ripk3                                                                          |
| <b>Storage</b>         | 2-8 °C for 6 months                                                            |
| <b>Note</b>            | For research use only                                                          |
| <b>Uniprot ID</b>      | <b>Q9QZL0</b>                                                                  |
| <b>Sensitivity</b>     | 0.094ng/ml                                                                     |
| <b>Expiration Date</b> | Please enquire.                                                                |



Mouse Ripk3  
(Receptor-interacting  
serine...