


Product Datasheet

Mouse EDA2R ELISA Kit (orb1292746)

Description	Mouse EDA2R (Ectodysplasin A2 Receptor)ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	EDA2R
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q8BX35
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
EDA2R(Ectodysplasin
A2 Receptor) E...