


Product Datasheet

Mouse Cxcr6 ELISA Kit (orb1292754)

Description	Mouse Cxcr6 (C-X-C chemokine receptor type 6) ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/ml
Target	Cxcr6
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9EQ16
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Mouse
Cxcr6 (C-X-
C
chemokine
receptor
ty...