


Product Datasheet

Mouse Cdnf ELISA Kit (orb1292755)

Description	Mouse Cdnf (Cerebral dopamine neurotrophic factor) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	Cdnf
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q8CC36
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse Cdnf
(Cerebral
dopamine
neurotroph...