


Product Datasheet

Mouse Hnrnpa2b1 ELISA Kit (orb1292761)

Description	Mouse Hnrnpa2b1 (Heterogeneous nuclear ribonucleoproteins A2/B1) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	Hnrnpa2b1
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
Hnrnpa2b1(Heterogeneous
nuclear ri...