


Product Datasheet

Mouse Prdx6 ELISA Kit (orb1292765)

Description	Mouse Prdx6 (Peroxiredoxin-6) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	Prdx6
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O08709
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
Prdx6(Peroxiredoxin-
6) ELISA Kit