


Product Datasheet

Mouse MX2 ELISA Kit (orb1292767)

Description	Mouse MX2 (Interferon-induced GTP-binding protein Mx2) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	MX2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9WVP9
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
MX2(Interferon-
induced GTP-
binding...