


Product Datasheet

Mouse MCPIP-1 ELISA Kit (orb1292774)

Description	Mouse MCPIP-1 (Monocyte Chemoattractant Protein-induced Protein 1) ELISA Kit
Reactivity	Mouse
Range	3.9-250pg/ml
Target	MCPIP-1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q5D1E7
Sensitivity	2.34pg/ml
Expiration Date	Please enquire.



Mouse MCPIP-1(Monocyte Chemoattractant P...