


Product Datasheet

Mouse Fundc1 ELISA Kit (orb1292780)

Description	Mouse Fundc1 (FUN14 domain-containing protein 1) ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/ml
Target	Fundc1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9DB70
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Mouse
Fundc1(FUN14
domain-
containing
pro...