


## Product Datasheet

### Mouse Efna4 ELISA Kit (orb1292781)

|                        |                                   |
|------------------------|-----------------------------------|
| <b>Description</b>     | Mouse Efna4 (Ephrin-A4) ELISA Kit |
| <b>Reactivity</b>      | Mouse                             |
| <b>Range</b>           | 31.25-2000pg/ml                   |
| <b>Target</b>          | Efna4                             |
| <b>Storage</b>         | 2-8 °C for 6 months               |
| <b>Note</b>            | For research use only             |
| <b>Uniprot ID</b>      | <b>O08542</b>                     |
| <b>Sensitivity</b>     | 18.75pg/ml                        |
| <b>Expiration Date</b> | Please enquire.                   |



Mouse  
Efna4(Ephrin-  
A4) ELISA Kit