


Product Datasheet

Mouse Enpp1 ELISA Kit (orb1292791)

Description	Mouse Enpp1 (Ectonucleotide pyrophosphatase/phosphodiesterase family member 1) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/ml
Target	Enpp1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P06802
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
Enpp1(Ectonucleotide
pyrophosphata...