


Product Datasheet

Mouse Arntl ELISA Kit (orb1292795)

Description	Mouse Arntl (Aryl hydrocarbon receptor nuclear translocator-like protein 1)ELISA kit
Reactivity	Mouse
Range	15.6-1000pg/ml
Target	Arntl
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9WTL8
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
Arntl(Aryl
hydrocarbon
receptor
nu...