


Product Datasheet

Mouse BCAT2 ELISA Kit (orb1292816)

Description	Mouse BCAT2 (Branched-chain-amino-acid aminotransferase, mitochondrial) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	BCAT2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O35855
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
BCAT2(Branched-
chain-amino-acid
am...