


Product Datasheet

Mouse BCAT1 ELISA Kit (orb1292817)

Description	Mouse BCAT1 (Branched-chain-amino-acid aminotransferase, cytosolic) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	BCAT1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P24288
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
BCAT1(Branched-
chain-amino-acid
am...