


Product Datasheet

Mouse TMPRSS2 ELISA Kit (orb1292822)

Description	Mouse TMPRSS2 (Transmembrane Protease, Serine 2) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	TMPRSS2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9JIQ8
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
TMPRSS2(Transmembrane
Protease, Se...