


Product Datasheet

Mouse Mcu ELISA Kit (orb1292831)

Description	Mouse Mcu (Calcium uniporter protein, mitochondrial) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	Mcu
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q3UMR5
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
Mcu(Calcium
uniporter
protein,
mit...