


Product Datasheet

Mouse ADAMTS8 ELISA Kit (orb1292834)

Description	Mouse ADAMTS8 (A disintegrin and metalloproteinase with thrombospondin motifs 8) ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	ADAMTS8
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P57110
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
ADAMTS8(A
disintegrin
and
metallop...