


Product Datasheet

Mouse Anti SARS-CoV-2 ELISA Kit (orb1292836)

Description	Mouse Anti SARS-CoV-2 (COVID) Spike RBD Alpha B.1.1.7 Variant (UK Variant) IgM ELISA Kit
Reactivity	Mouse
Range	0.781-50ng/ml
Target	Anti SARS-CoV-2 UK Variant IgM
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	0.469ng/ml
Expiration Date	Please enquire.



Mouse Anti SARS-CoV-2 (COVID) Spike RBD ...