


Product Datasheet

Mouse ALDOB ELISA Kit (orb1292838)

Description	Mouse ALDOB (Fructose-bisphosphate aldolase B)ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	ALDOB
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q91Y97
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
ALDOB(Fructose-
bisphosphate
aldola...