


Product Datasheet

Mouse Sfrp4 ELISA Kit (orb1292854)

Description	Mouse Sfrp4 (Secreted frizzled-related protein 4) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	Sfrp4
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9Z1N6
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
Sfrp4(Secreted
frizzled-
related pr...