


Product Datasheet

Mouse UCP-1 ELISA Kit (orb1292857)

Description	Mouse UCP-1 (Mitochondrial brown fat uncoupling protein 1) ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/ml
Target	UCP-1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P12242
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Mouse UCP-1 (Mitochondrial brown fat uncoupling protein 1) ELISA Kit