


## Product Datasheet

### Human MBD2 ELISA Kit (orb1292864)

|                        |                                                               |
|------------------------|---------------------------------------------------------------|
| <b>Description</b>     | Human MBD2 (Methyl-CpG-binding domain protein 2)<br>ELISA Kit |
| <b>Reactivity</b>      | Human                                                         |
| <b>Range</b>           | 15.625-1000pg/ml                                              |
| <b>Target</b>          | MBD2                                                          |
| <b>Storage</b>         | 2-8 °C for 6 months                                           |
| <b>Note</b>            | For research use only                                         |
| <b>Sensitivity</b>     | 9.375pg/ml                                                    |
| <b>Expiration Date</b> | Please enquire.                                               |



Human  
MBD2(Methyl-  
CpG-binding  
domain pro...