


Product Datasheet

Human MAPKAPK2 ELISA Kit (orb1292866)

Description	Human MAPKAPK2 (MAP kinase-activated protein kinase 2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	MAPKAPK2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P49137
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
MAPKAPK2(MAP
kinase-
activated
prot...