


Product Datasheet

Human CACNA1C ELISA Kit (orb1292872)

Description	Human CACNA1C (Voltage-dependent L-type calcium channel subunit alpha-1C)ELISA Kit
Reactivity	Human
Range	78.125-5000pg/ml
Target	CACNA1C
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q13936
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Human
CACNA1C
(Voltage-
dependent
L-type ...