


Product Datasheet

Human HSD3B7 ELISA Kit (orb1292899)

Description	Human HSD3B7 (3 beta-hydroxysteroid dehydrogenase type 7) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	HSD3B7
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9H2F3
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
HSD3B7(3
beta-
hydroxysteroid
dehyd...