


Product Datasheet

Human GOLPH3 ELISA Kit (orb1292917)

Description	Human GOLPH3 (Golgi phosphoprotein 3) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	GOLPH3
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9H4A6
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
GOLPH3(Golgi
phosphoprotein
3) ELI...