


Product Datasheet

Human GOLGA2 ELISA Kit (orb1292918)

Description	Human GOLGA2/GM130 (Golgin subfamily A member 2) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	GOLGA2/GM130
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q08379
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
GOLGA2/GM130(Golgin
subfamily A me...