


Product Datasheet

Human GAR1 ELISA Kit (orb1292925)

Description	Human GAR1 (H/ACA ribonucleoprotein complex subunit 1) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	GAR1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9NY12
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
GAR1(H/ACA
ribonucleoprotein
compl...