


Product Datasheet

Human GABRA1 ELISA Kit (orb1292930)

Description	Human GABRA1 (Gamma-aminobutyric acid receptor subunit alpha-1) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	GABRA1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P14867
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
GABRA1(Gamma-
aminobutyric
acid rec...