


Product Datasheet

Human CHRNA7 ELISA Kit (orb1292968)

Description	Human CHRNA7 (Neuronal acetylcholine receptor subunit alpha-7) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/ml
Target	CHRNA7
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P36544
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
CHRNA7(Neuronal
acetylcholine
rece...