


Product Datasheet

Human anti-SARS-CoV2 ELISA Kit (orb1293005)

Description	Human anti-SARS-CoV2 (S-RBD) (Omicron,B.1.1.529) IgG ELISA Kit
Reactivity	Human
Range	3.9-250ng/ml
Target	anti-SARS-CoV2(S-RBD) (Omicron,B.1.1.529)
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	2.34ng/ml
Expiration Date	Please enquire.



Human anti-SARS-CoV2(S-RBD) (Omicron, B....