


Product Datasheet

Human CHRNA2 ELISA Kit (orb1293007)

Description	Human CHRNA2 (Neuronal acetylcholine receptor subunit alpha-2) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	CHRNA2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q15822
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
CHRNA2(Neuronal
acetylcholine
rece...