


Product Datasheet

Human Anti-VCL ELISA Kit (orb1293030)

Description	Human Anti-VCL (Anti-Vinculin Antibody)ELISA Kit
Reactivity	Human
Range	1.563-100ng/ml
Target	Anti-VCL
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	HRP□1□100
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Human
Anti-
VCL(Anti-
Vinculin
Antibody)
E...