


## Product Datasheet

### Human FBLN7 ELISA Kit (orb1293034)

|                        |                                   |
|------------------------|-----------------------------------|
| <b>Description</b>     | Human FBLN7 (Fibulin-7) ELISA Kit |
| <b>Reactivity</b>      | Human                             |
| <b>Range</b>           | 31.25-2000pg/ml                   |
| <b>Target</b>          | FBLN7                             |
| <b>Storage</b>         | 2-8 °C for 6 months               |
| <b>Note</b>            | For research use only             |
| <b>Uniprot ID</b>      | <b>Q53RD9</b>                     |
| <b>Sensitivity</b>     | 18.75pg/ml                        |
| <b>Expiration Date</b> | Please enquire.                   |



Human  
FBLN7(Fibulin-  
7) ELISA Kit