


Product Datasheet

Human DYRK1A ELISA Kit (orb1293035)

Description	Human DYRK1A (Dual specificity tyrosine-phosphorylation-regulated kinase 1A) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	DYRK1A
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q13627
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
DYRK1A(Dual
specificity
tyrosine-p...