


Product Datasheet

Human CHD8 ELISA Kit (orb1293036)

Description	Human CHD8 (Chromodomain-helicase-DNA-binding protein 8) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	CHD8
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9HCK8
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
CHD8(Chromodomain-
helicase-DNA-bin...