


Product Datasheet

Human CEFIP ELISA Kit (orb1293042)

Description	Human CEFIP (Cardiac-enriched FHL2-interacting protein) ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	CEFIP
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q711Q0
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
CEFIP(Cardiac-
enriched
FHL2-intera...