


Product Datasheet

Human ZNF207 ELISA Kit (orb1293051)

Description	Human ZNF207 (BUB3-interacting and GLEBS motif-containing protein ZNF207) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/ml
Target	ZNF207
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O43670
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
ZNF207(BUB3-
interacting
and GLEBS ...