


Product Datasheet

Human GRIN2D ELISA Kit (orb1293064)

Description	Human GRIN2D (Glutamate receptor ionotropic, NMDA 2D)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	GRIN2D
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O15399
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
GRIN2D(Glutamate
receptor
ionotrop...