


Product Datasheet

Human VISTA ELISA Kit (orb1293070)

Description	Human VISTA (V-type immunoglobulin domain-containing suppressor of T-cell activation) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	VISTA
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9H7M9
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human VISTA
(V-type
immunoglobulin
domai...