


Product Datasheet

Human FGD4 ELISA Kit (orb1293076)

Description	Human FGD4 (FYVE, RhoGEF and PH domain-containing protein 4)ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	FGD4
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q96M96
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
FGD4(FYVE,
RhoGEF and
PH domain-
CO...