


Product Datasheet

Human PRKAR2A ELISA Kit (orb1293077)

Description	Human PRKAR2A (cAMP-dependent protein kinase type II-alpha regulatory subunit) ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	PRKAR2A
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P13861
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
PRKAR2A
(cAMP-
dependent
protein ki...